



PRODUCTS


Different sizes of DR-MH Series

Ultrasonic Cleaner P Series	Model	DR-MH13	DR-MH20	DR-MH22	DR-MH30	DR-MH40
Ultrasonic Cleaner DS Series	Tank Size (L)	1.3	2	2.2	3	4
Ultrasonic Cleaner LD Series	Tank Internal Dimension (mm)	150×138×65	150×138×100	238×138×65	238×138×100	302×152×100
Ultrasonic Cleaner LQ Series	Operation Frequency (KHz)	40	40	40	40	40
Ultrasonic Cleaner MH Series	Ultrasonic Power	60	60	120	120	120
Ultrasonic Cleaner M&MS Series	Heating Power	60	100	100	150	220
Derui Ultrasonic Cleaner-DSP S	Drainage Valve	No	No	Optional	Optional	Optional
Ultrasonic Homogenizer	Carrying Handle	No	No	Yes	Yes	Yes
Ultra Transducer	Cover	Yes	Yes	Yes	Yes	Yes
Accessories	Unit Weight (Kg)	3.5	3.5	4	4.5	5.8
	Unit /Carton	4	4	4	4	2

Carton size(cm)	53x53x36	53x53x36	52x41x67	52x38x67	60x45x38
Carton Weight(kg)	13.8	12.6	15.8	18	11.6

Model	DR-MH100	DR-MH120	DR-MH130	DR-MH150	DR-MH220
Tank Size (L)	10	12	13	15	22
Tank Internal Dimension (mm)	302×239×150	490×138×150	302x239x200	325x295x150	502×298×150
Operation Frequency (KHz)	40	40	40	40	40
Ultrasonic Power	240	300	300	360	360
Heating Power	310	380	420	470	470
Drainage Valve	Yes	Yes	Yes	Yes	Yes
Carrying Handle	Yes	Yes	Yes	Yes	Yes
Cover	Yes	Yes	Yes	Yes	Yes
Unit Weight (Kg)	6.5	6.8	7.2	8.6	10.1
Unit /Carton	1	1	1	1	1
Carton size(cm)	47x36x43	67x26x47	48x37x45	48x43x47	67x43x42
Carton Weight(kg)	8.5	8	9.5	12	13.5

Feature Illustration



MH Series Feature

- The stainless steel tank used and the reliable technology prolong the service life of the units
- With mechanical timer & heater control
- Ultrasonic Timer 0~60 Minute, infinity
- Temperature adjust from 25°C to 80°C
- High-performance transducer system

Here is a partial list of some of the things we have helped clean



- ♦ Jewellery: gold, silver, platinum and fancy jewellery, Necklaces, rings, Bracelets etc.
 - ♦ Watch and glasses: watch cases, mechanical parts and metal bracelets, waterproof wristwatches, watch chains, gl etc.
 - ♦ Dental laboratory: instruments, removal of cement and plaster; cleaning of dental mould spoons, drills and pros crowns and other prosthetic materials.
 - ♦ CD DVD: Any CD VCD DVD pre-recorded of CD-R/RW
 - ♦ Office equipment: Printer heads, Fountain pen nibs Wax or Rubber stamps etc.
 - ♦ Mental Items: Metal cutlery (knives/forks etc), old coins, mental badges, valves, machine nozzles, small metal pa
 - ♦ Laboratory and Medical sector: laboratory instruments made of glass, plastic or metal, surgical, medical and of i instruments without damaging the material, or cleaning of rigid endoscopes and flexible endoscopic accessories
- © 2012 Derui International Limited. All rights reserved