



HY-8 Large Capacity Multi-use shaker

PRODUCT DETAIL

Introduction:

The shaker is one with more advanced technique and reliable quality in its like. It has advantages of stable operation, low noise, stepless adjustable shaking speed and wide adjustable range. The spring multi-purpose tray can match with diversified flasks to meet needs from different customers. It is an optimal equipment in research and production in biochemical industry, biological engineering, teaching and studying, microbiology and medical science and so on.

Technical Data:

Model	Shaking Frequency	Load	Shaking Mode	Power Supply	Overall Size
HY-1/HY-3	Start-300	20KG	Vertical/Reciprocating	220V 50HZ	50X30X34CM
KS/HY-2A/HY-4	Start-300	30KG	Reciprocating	220V 50HZ	40X26X30CM
HZ-200/HZ-300 HY-B/HY-5	Start-300	30KG	Rotary	220V 50HZ	40X40X24CM
HY-6	Start-250	60KG	Reciprocating/Rotary	220V 50HZ	62X47X58CM
HY-8	Start-260	100KG	Reciprocating/Rotary	220V 50HZ	76X62X32CM
TY-80R/S	0-120	5KG	Rotary/Swing	220V 50HZ	24X32X16CM