

Petri Dish Plastic



Produced with high quality polystyrene (PS), Petri Dishes have smooth surface and optimal transparent body. With EO sterilization, apply for bacteria & cell culture. Other-more the stackable design makes storage easy.

- To prevent condensation by the three venting (Triple)
- Sterilized by Ethylene Oxide (ETO) density 1 kg/m³
- Clear color, non-autoclave
- Maximum temperature: 85°C (185°F)
- Brittleness temperature: -20°C (-4°F)

CODE.	Description	SIZE	Weight	Packaging
PDB9015	Sterilized Petri Dish Plastic, Polystyrene (Premium)	90x15mm	17.5 G	Plastic/Paper pouch bag 500 Pairs (10Pairs/bag)

Polystyrene has high strength and impact resistance and excellent compatibility with weak acids, alcohol and bases.