

MICRO QC

Plastic Petri Dish Sterillized

Produced with high quality polystyrene (PS), Petri Dishes have smooth surface and optimal transparent body. With EO sterilization, apply for bacteria & cell culture. Othermore the stackable design makes storage easy.

- Maximum temperature: 80°C (185°F)
- Brittleness temperature: -20°C (-4°F)
- To prevent condensation by the three venting (Triple)
- Sterilized by Ethylene Oxide (ETO) density 1 kg/m³
- Clear color, non-autoclave

CODE.	Description	Weight	SIZE	Packaging
PDMQ-9015	Sterilized Petri Dish Plastic, Polystyrene (Premium)	14G	90x15mm	Plastic/Paper pouch bag 500Pairs (10 Pairs/bag)



⊘ STERILE UNLESS PACKAGE IS OPENED OR DAMAGED