



Personalized PDF Catalog
Catalogue Generated 21 January, 2019

Four-way racks for microtubes and test tubes

Supplier: VWR Collection



PP.

- »Side 1: 4×30 mm diameter holes for 50 ml centrifuge tubes
- »Side 2: 6×17 mm diameter holes for 15 ml centrifuge tubes
- »Side 3: 32×12 mm diameter holes for 1,5/2,0 ml microcentrifuge tubes
- »Side 4: 32×8 mm diameter holes for 0,5 ml microcentrifuge tubes

With four sides that can be loaded independently. Create as large a rack as needed by linking together multiple racks. Practical and compact. Very good chemical resistance. High thermal resistance, autoclavable.

ORDER

VWR Catalog Number	Quantity	Colour	W×D×H	Pk
VWRI82010-750		Assorted colours	175×95×51 mm	5

Supplier Article Number:: 82010-750

VWRI525-0952		Assorted colours	174×95×52 mm	1 SET
VWRI525-0951		Blue	174×95×52 mm	1
VWRI211-0201		Green	174×95×52 mm	1
VWRI82024-460		Green	175×95×51 mm	5

Supplier Article Number:: 82024-460

VWRI82010-752		Natural	175×95×51 mm	5
---------------	--	---------	--------------	---

Supplier Article Number:: 82010-752

VWRI525-0950		Natural	174×95×52 mm	1
VWRI211-0200		Orange	174×95×52 mm	1
VWRI211-0202		Pink	174×95×52 mm	1
VWRI211-0203		Yellow	174×95×52 mm	1

This product is no longer available. Alternatives may be available by searching with the VWR Catalog Number listed above. If you need further assistance, please contact us at +65 6505 0760.