

Biochemistry Incubator



BJPX-BI

BJPX-BII

BJPX-BIII

- Features:
- * BJPX-BI series: LED display.
 - * Polished stainless steel chamber.
 - * CFC-free eco-friendly refrigerant.
 - * BOD socket integrated inside for additional small equipment.
- * BJPX-BII & BIII series: LCD display
 - * Round angle structure, easy to clean.
 - * Microprocessor PID temperature control, accurate and reliable.
- * Built-in printer for BJPX-BIII series
 - * Universal castors for easy movement.

- Safety Functions:
- * Temperature protection threshold value can be set.
 - * Auto-off protection and alarm for abnormal situation.
 - * Parameters can be saved automatically when power off or system error.

- Optional Accessories:
- * Test port.

Technical Parameters:

Model	BJPX-B80I	BJPX-B150I	BJPX-B200I	BJPX-B250I	BJPX-B300I	BJPX-B400I
	BJPX-B80II	BJPX-B150II/III	BJPX-B200II/III	BJPX-B250II/III	BJPX-B300II/III	BJPX-B400II/III
Capacity	81L	160L	202L	251L	304L	397L
Temp. Range	5~50°C (type II/III 0~65°C)					
Temp. Fluctuation	±1°C (type II & III: heating: ±0.5°C to ±1°C)					
Temp. Uniformity	±2°C					
Shelves	2 pcs/adjustable	3 pcs/adjustable				
Time Range	0~9999min					
Refrigerant	R134a (CFC free, eco-friendly)					
BOD socket	1pc					
Consumption	650W	710W	760W	810W	960W	1280W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)					
Internal Size(mm)	400*370*550	450*420*850	450*450*1000	480*490*1070	520*500*1170	580*540*1270
External Size(mm)	540*580*1075	598*633*1471	590*659*1525	620*690*1595	660*702*1695	720*747*1795
Package Size(mm)	700*710*1260	760*810*1660	790*740*1700	830*770*1750	790*780*1850	850*830*1960
Gross Weight(kg)	87	116	146	147	157	170

Biochemistry Incubator

- Features:
- * LCD display.
 - * Polished stainless steel chamber.
 - * Round angle structure, easy to clean.
 - * CFC-free environmental friendly refrigerant.
 - * Double-layer glass observation window design.
 - * Microprocessor PID temperature control: accurate and reliable.
 - * Universal castors for easy movement (except the model BJPX-B70/100GK)



LCD display



BJPX-B70/100GK



BJPX-B150/250/350/450GK

- Safety Functions:
- * Temperature protection threshold value can be set.
 - * Function of self-diagnosis and alarm for abnormal situation.
 - * Parameters can be saved automatically when power off or system error.

- Optional Accessories:
- * Test port.
 - * BOD socket.
 - * RS485 interface or built-in printer.

Technical Parameters:

Model	BJPX-B70GK	BJPX-B100GK	BJPX-B150GK	BJPX-B250GK	BJPX-B350GK	BJPX-B450GK
Capacity	70L	100L	150L	250L	350L	450L
Temp. Range	4~60℃					
Temp. Fluctuation	±0.5℃					
Temp. Precision	0.1℃					
Ambient Temp.	+4~35℃					
Shelves	2 pcs/ adjustable		3pcs/ adjustable			4pcs/ adjustable
Time Range	1~9999min					
Refrigerant	R134a (CFC free, eco-friendly)					
Working Type	Continuous					
Consumption	660W	680W	690W	730W	780W	1320W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)					
Internal Size(mm)	450*350*450	510*380*510	505*380*800	550*450*1050	550*550*1150	700*550*1200
External Size(mm)	550*550*890	610*580*940	630*600*1450	680*660*1710	680*770*1870	830*750*1860
Packing Size(mm)	680*680*1050	740*710*1140	760*710*1610	790*770*1860	810*900*2020	880*950*2020
Gross Weight(kg)	60	100	120	150	160	180