

Acid & alkali storage cabinet

PRODUCT INTRODUCTION

SAFOO series high performance safety cabinet strictly corresponds to OSHA29 CFR1910.106 and NFPA CODE 30 STANDARD, at the same time, it passes third party EU CE authorized certification.

Strong acid&alkali storage cabinet-product feature



Cabinet door: apply 15mm PP(polypropylene) material.

PP laminate ,Separate extraction, Easy to clean

Bridge doorknob: PP(polypropylene) material with anti-corrosive performance(Color selectable: azure blue, Porcelain white, Light grey ,Dark grey)



In order to provide safety, double-lock structure can enhance the safety factor.



Cabinet body apply 8mm PP(polypropylene) plate with Seamless welding, beautiful appearance, with more than 10 years of service life.



Warning labels
Doorknob is marked with a marked corrosive warning.



Hing: apply injection molding PP material ,with good anti-corrosive performance. Color is selectable.



Apply PP(polypropylene) material, Longer service life, Quality guaranteed.

Acid base storage cabinet

Model	Parameter	Configure
SF004P	Dimension(mm): 560 × 430 × 430mm Specification: 4 gallon Standard laminate: 1	Single door ,double-lock
SF012P	Dimension(mm): 890*590*460mm Specification: 12 gallon Standard laminate: 1	Single door ,double-lock
SF028P	Dimension(mm): 900*910*600mm Specification: 28 gallon Standard laminate: 1	Double door,double-lock
SF030P	Dimension(mm): 1120*1090*460mm Specification: 30 gallon Standard laminate: 1	Double door,double-lock
SF045P	Dimension(mm): 1650*1090*460mm Specification: 45 gallon Standard laminate: 2	Double door,double-lock
SF060P	Dimension(mm): 1650*860*860mm Specification: 60 gallon Standard laminate: 2	Double door,double-lock
SF090P	Dimension(mm): 1650*1090*860mm Specification: 90 gallon Standard laminate: 2	Double door,double-lock