

The high-performance

## CHEMICAL REAGENTS STORAGE CABINET

with its advanced features and great ergonomic design



### SWN series models

With the advanced Magnetic Gasket on the door frame, it enables every user to easily open and close the door as well as excellent effect for sealing.

One-touch replacement system for filter gives very good merit, easiness in maintenance.

User-oriented safety design and high-effectiveness powerful motor with low noise

Anti-corrosive construction for the storage room made of ABS injection molding and it features excellent protection and high-effectiveness against corrosi.



GENERAL SPECIFICATIONS		SWN 2001 / 2002
Material	Interior	Steel Plate with Powder Coating
	Exterior	ABS Injection Molding
Filter		SHIN SAENG FILTER(Prefilter+Carbon Filter+Prefilter)
Fan		High Efficiency and Silent Fan(Twin Fan System)
Noise Level		Below 60 decibels
Controller		LCD Display(Digital)
Lamp		L.E.D(High Luminance Lamp) × 2ea
Handle		Moving Lock Handle
Caster		Packing Caster incuding function for stop × 4ea
Supply Power&Frequency		AC230V, 50/60Hz

GENERAL SPECIFICATIONS	SWN 2001	SWN 2002
External Dimension(L×D×H)	1250×760×2005(mm)	900×760×2005(mm)
Shelves	12 ea	6 ea

Our ceaseless R&D and improvements may result in specification changes at any time.  
Different size, specification and color is always available in accordance with customer's request.

For details on these superior products or any other inquiries, please feel free to contact your local sales representative or you can contact us by e-mail at [contact@shinsaeng.com](mailto:contact@shinsaeng.com)