

BRAND: PULISITE @ **HUAZHI**

MANUFACTURE: HUAZHI (FUJIAN) ELECTRONIC TECHNOLOGY CO.,LTD

Focusing on the development and manufacture of excellent balances  
and laboratory instruments for years



## PTT-A Precision Scale (0.01g)

### Application Mode

- ✧ Percentage Weighing
- ✧ Selectable weighing and count mode

**Optional with RS232:** Weighing functions are selectable on PC. Measurements can be input into EXCEL, WORD and other files on PC.

### Characteristic

- ◆ Loadcell Sensor
- ◆ External Calibration
- ◆ Six level anti-vibration adjustable
- ◆ Weighing stability time adjustable
- ◆ More than 14 unit shield: g, kg, ct...
- ◆ Default weighing mode automatically
- ◆ Optional with rechargeable battery
- ◆ Optional with RS232 to communicate with external devices, like PC and printer
- ◆ Optional with rechargeable battery
- ◆ Optional with external display
- ◆ Optional with underhook


### External CAL (With CAL weights) Loadcell Sensor

PTT is a standard product with the function of AC-DC which is convenient for outdoor users. It has LCD display and is optional with external display. It is the most cost-effective balance.


Optional with  
windshield



### Characteristic



### Application Mode



Model	Capacity(g)	Readability	Repeatability	Linearity(g)	Pan(mm)	Response Time	Housing(mm)
PTT-A3200	3200g	0.01g	±0.02g	±0.03g	156×156	2-3 sec.	295×208×88