


Levo Pipette Controller

Features

- Reduce hands-on time
- Compatible with glass or plastic serological pipettes (0.1 – 100mL)
- 3µm replaceable hydrophobic filter
- Easy to clean and maintain



Mini Pipette Features

- Ergonomic design provides excellent operation experience;
 - Smart body size ideal to be packed into IVD assay kits;
- Compatible with universal 200µL/300µL tips;

Volumes 300µL 200µL 100µL 80µL 75µL 50µL 25µL OEM

Levo E Pipette Pump

Designed for single-handed operation, Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The feature of side lever enables complete liquid dispensing, without a droplet remaining.



Features

- Volume capacity 2mL, 10mL and 25mL
- Colour-coded by volume with green, blue and red
- Thumb wheel guarantees precise operation
- Resistant to acids and alkalis
- Easy for cleaning and maintenance