



Levo Plus Pipette Filler

Features

- LCD display to observe battery level, aspiration/ dispensing speed.
- Light in weight and ergonomic design.
- Compatible with volumetric and serological pipettes.
- Compatible pipette volume range of 0.1-100mL.
- Integrated with replaceable, hydrophobic filter of 0.45µm.
- Different speeds of selection for aspiration and dispensing.
- Lithium-ion battery enables more than 8hours of intermittent use.



Specifications

Aspiration Speed	25mL<5s (8 shift)
Dispensing Speed	Motor (8 shift)/Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45µm Hydrophobic
Weight	200g