

VERTICAL PRESSURE STEAM STERILIZER

(Digital Display Automation)



Characteristics:

- Fully stainless high-alloy steel structure
- Digital display of working status, touch of key
- Auto discharge the cool air, and steam discharging automatically after sterilization
- Automatically shut off with beep reminding after sterilization
- With two stainless steel sterilizing baskets
- Easy to operate, safe and reliable

Technical Data:

Model technical data	LS-35LD(automatic)	LS-50LD(automatic)	LS-75LD(automatic)	LS-100LD(automatic)
Chamber volume	35L(φ318×450)mm	50L(φ340×550)mm	75L(φ400×600)mm	100L(φ440×650)mm
working pressure	0.22MPa			
Working temperature	134℃			
Max working pressure	0.23 Mpa			
Heat average	±1℃			
Timer	0~99min or 0~99hour59min			
Adjustment of temperature	105~134℃			
Power	2.5Kw/AC220V.50Hz	3Kw /AC220V.50Hz	4.5KW AC220V.50Hz	
Overall dimension	470×450×840(mm)	520×500×970(mm)	550×530×990(mm)	580×560×1060(mm)
Transport dimension	540×540×970(mm)	580×580×1090(mm)	610×610×1110(mm)	640×640×1170(mm)
G.W/N.W	61Kg/ 46Kg	73Kg/56Kg	87Kg/67Kg	95Kg/73Kg