



Uniclave Series Autoclave

Type A: Fully automatic, basic Model

Type R: Drying function based on Type A

◆ FunctionTable

Model	FD36A	FD50A	FD80A	FD100A
	FD36R	FD50R	FD80R	FD100R
Capacity(L)	36	50	80	100
Outer dimension(mm)	460x542x960	460x542x1070	620x680x1042	620x680x1222
Sterilization chamber size(mm)	325x460	325x670	400x709	400x889
Rated power supply	220V 50Hz/60Hz			
Rated power(W)	2300	2600	4600	4600
Working environment	5℃-40℃, relative humidity 10%-85%			
Chamber material	SUS 304 stainless steel			
Sterilization temperature	105-135℃			
Sterilization time	1-300 min			
Drying Time	1-300 min(only available for Type R)			
Range of pressure gauge	0-0.4Mpa			
Controller	Microcomputer control system			
Mode	Solid mode, liquid mode, drying function of Type R			
Safety Device	Self-induction interlocking device, over temperature protection system, dry scorch protection, over pressure protection, safety valve, over current and short circuit protection, cover checking system, electric leakage protection, scald protection, automatic failure-inspection system			
Standard accessories	Stainless steel bakets, waste bottle, water plate, drainage hose			
Optional parts	Printer, printing set, stainless steel pails			

