


- Laboratory Shaker/Flask Shaker
- 1. Bi-layer, huge capacity, 2 shake mode
 - 2. ISO13485/14001/9001, SGS audited factory
 - 3. Since 1984

The shaker is one with more advanced technique and reliable quality in its like. It has advantages of stable operation, low noise, stepless adjustable shaking speed and wide adjustable range. The spring multi-purpose jigger can match with diversified flasks to meet needs from different customers. It is an optimal equipment in research and production in biochemical industry, biological engineering, teaching and studying, microbiology and medical science and so on.



Specifications of HY-6 Double Stage Laboratory Shaker/Flask Shaker	
Power Supply	AC 220V/110V 50Hz/60Hz
shaking amplitude	20mm(convoluting), 40mm(reciprocating)
Working method	convoluting or reciprocating customizing
Timer	0-120min
Working load	60kg
shaking frequency	0~250times/min stepless regulation (digital display velocity)
Capacity	500ml*20*2stage or 250ml*25*2stage customizing
Oscillation tray	Bilayer 62X47mm
Overall size	62X47X58cm