

Mini Centrifuge

Mini-15K High Speed Centrifuge

■ Features:

1. Extremely quiet operation
2. Easy-to-follow LCD display
3. Capacity: 12 x 1.5/2.0ml tubes
4. Separate short-spin key
5. Low sample heating (only 12°C after 30min, at max.speed)
6. Autoclavable rotor (121°C, 20min)
7. RMP/RCF setting as required
8. Includes standard rotor made of anodized aluminum
9. Timer can be set up to 99min
10. The small size of the unit makes it convenient for benchtop use or in a hood or in refrigerated environments

■ Specification

Model	Mini-15K
Power requirement	100W
Relative centrifugal force (RCF)	14,000 x g
Max.rotational speed	14,500rpm
Timer	20sec~99min
Acceleration time to max. speed	15s
Braking time from max. speed	15s
Dimensions (DxWxH)	270x230x155mm



Code	Description
AS-08110-00	Mini-15K AC100V or 240V 50~60Hz 3A
AS-08041-01	Adapters for 0.5ml tubes
AS-08041-02	Adapters for 0.2ml tubes

Mini-6KC Centrifuge

Mini-6KC Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm & 6000rpm). The electronic brake can reach quick deceleration to reduce handling time.

■ Features:

1. New rotor design for increased sample capacity
2. Quick release rotor hub allows for easy rotor exchange
3. The electronic braking system for rapid sample processing
4. Two speed options:4,000 & 6,000rpm

■ Specification

Model	Mini-6KC
Capacity	8x1.5/2.0 ml tube rotor & 4 place PCR Strip rotor
Two Speed options	4,000rpm & 6,000rpm
Maximum RCF	1,200 x g (4,000rpm) / 2,000 x g (6,000rpm)
Dimensions	135 x 147 x 107mm (D x W x H)
Weight	0.7kg
Electrical	AC100~240V, 50/60Hz
Code	Description
AS-08020-00	Mini-6KC AC100V~240V 50~60Hz 1A
AS-08021-01	8-1.5/2.0ml angle rotor
AS-08021-02	4 - 8x0.2ml strip rotor
AS-08061-05	Adapters for 0.5~0.6ml tubes
AS-08061-06	Adapters for 0.2ml tubes



Two speed options



Quick release rotor