

Multi-use shaker

The shaker is one with more advanced technique and reliable quality in its like. It has advantages of stable operation, low noise, stepless adjustable shaking speed and wide adjustable range. The spring multi-purpose tray can match with diversified flasks to meet needs from different customers. It is an optimal equipment in research and production in biochemical industry, biological engineering, teaching and studying, microbiology and medical science and so on.



HY-4(KS)
Speed - adjustable Multipurpose Shaker



TY-80R/S
Decolorization Shaker



HY-2.4(A)
Speed - adjustable Multipurpose Shaker



MM-2
MM-2 Micro Shaker



HZ-300、HZ-200
HZ-300、HZ-200 Rotary Shaker



HY-B2、HY-5
HY-B2、HY-5 Rotary Shaker



HY-3
HY-3 Multi - function Shaker



HY-6
HY-6 Bi - deck Shaker



HY-8
HY-8 Large capacity Shaker

Technical Data

Model	Name	Shaking Frequency (rpm)	Load	Shaking Mode	Power Supply	Overall Sizes
HY-1、HY-3	multi-function shaker	start-300	20KG	↑↓ —	220V、50HZ	50×30×34cm
KS、HY-2A、HY-4	speed adjustable multipurpose shaker	start-300	30KG	—	220V、50HZ	40×26×30cm
HZ-200、HZ-300 HY-B、HY-5	rotary shaker	start-300	30KG	↻	220V、50HZ	40×40×24cm
HY-6	bi-layer shaker	start-250	60KG	— ↻	220V、50HZ	62×47×58cm
HY-8	large size speed adjustable shaker	start-260	100KG	— ↻	220V、50HZ	76×62×32cm
TY-80R/S	decolorization shaker	0-120	5KG	↻	220V、50HZ	24×32×16cm