

Digital Thermostatic Water Bath

High quality stainless steel material, heat resisting and anticorrosion, operational temperature controller and bi-integral A/D converter, digital display, durable.



HH-S₁、HH-S₂
HH-S₁、HH-S₂ Thermostatic water bath



HH-S₄
HH-S₄ Thermostatic water bath




HH-S₆
HH-S₆ Thermostatic water bath



HH-S₈
HH-S₈ Thermostatic water bath



HH-S₁₁、4、6、8
Single – row Thermostatic Water Bath



HH-S₁₀、12、15
HH-S₁₀、12、15 Large-capacity water bath

Technical Data

Model	HH-S ₁	HH-S ₂	HH-S _{11.4}	HH-S _{11.6}	HH-S ₄	HH-S ₆	HH-S ₈	HH-W420	HH-W600	DK-8D	W201B
Style	single-row	single-row	single-row	single-row	double-row	double-row	double-row	single-row	double-row	2. 3. 4row	single bath
Heating Power	300W	500W	1200W	1500W	1200W	1500W	1800W	800W	2000W	300W×X	800W
Temperature Control	RT—100℃										
Fluctuation of WaterTemperature	≤0.5℃										
Temperature Accuracy	≤0.5℃										
Working room(cm ²)	14×15×15	30×15×15	60×18×12	90×18×12	32×30×15	50×30×15	62×32×12	42×18×18	60×30×18	14×12×12	Φ 26×16