

NEW

# MS-T-S15

## 15-Position magnetic stirrer

### Features

- Efficient, simple, safe, multi-position agitator is the ideal choice for sample testing.
- Simple and easy operation panel: 15 stirring positions are controlled by 3 independent channels, all channels are available at different speeds and times.
- Brushless DC motor and powerful magnetic ensure fast and accurate stirring speed.
- Aluminum laser engraving disk surface, skid and corrosion proof anti-abrasive.



Specifications	MS-T-S15
Work plate Dimension [WXD]	460x304mm
Work plate material	Aluminum
Motor type	DC brushless
Stirring positions	15
Max. stirring quantity [HzO]	0.4Lx15
Max. magnetic bar[length]	40mm
Speed control accuracy	± 10rpm
Speed 、 Time display	LCD
Speed range	0–1500rpm
Timer function	0min–99h59min
Protection class	IP21
Power	40w
Voltage ,Frequency	110V/220V50Hz/60Hz
Dimension [DxWxH] [mm]	615x307x90mm
Weight	9.2kg
Permissible ambient temperature and humidity	5–40℃, 80%RH



0.4Lx15



0–1500rpm



Timer



Skid resistance



Multiple position



Brushless DC motor