

Microcentrifuge Tubes

Microcentrifuge Tubes

Polypropylene (PP) or translucent Polyethylene (LD-PE).

- autoclavable with opened cap at 121 °C, 20 min, (except 5615400)
- natural color or different colors for comfortable identification
- tight sealing cap
- safe and easy to be opened and closed
- 5615005 and 5615106: free of DNase-, RNase-, pyrogens-, Human DNA-, PCR inhibitors



CE IVD*

Volume ml	Fig. No.	Material	Cap	Color	Packaging	Quantity per Pack	Order No.
0.4 Beckman	1	PE	■	natural	1 x 1000	1000	56 15 400
0.5	2	PP	■	natural	1 x 1000	1000	56 15 030
1.5	3	PP	■	natural	2 x 500	1000	56 15 000
1.5	—	PP	■	yellow	2 x 500	1000	56 15 001
1.5	—	PP	■	green	2 x 500	1000	56 15 002
1.5	—	PP	■	red	2 x 500	1000	56 15 003
1.5	—	PP	■	blue	2 x 500	1000	56 15 004
1.5	3	PP	■	natural	1 x 1000	1000	New 56 15 005
1.5	4	PP	—	natural	1 x 1000	1000	56 15 020
2.0	5	PP	■	natural	1 x 1000	1000	56 15 101
2.0	5	PP	■	natural	1 x 1000	1000	New 56 15 106

*except 5615030

Safety-Cap Microcentrifuge tubes

Polypropylene (PP).

The new cap sealing of test tubes ensures an extreme tightness which is essential for applications with valuable samples, radioactive materials and DNA.

- autoclavable with opened cap at 121 °C, 20 min
- 0,5 ml and 1,5 ml centrifugable up to 30.000g
- 2,0 ml centrifugable up to 20.000g
- free of DNase-, RNase-, PCR Inhibitors
- free of endotoxin-, human DNA and ATP
- produced in microbiologically controlled clean room ISO Class 8 146441-1
- optimal tightness due to tight serration at the cap and tight channel at the tube; thereby a low evaporation is guaranteed
- secured shut due to form-closed cap lock; accidental opening of the lid during incubation is not possible
- high chemical resistance
- labeling areas on the cap and tube
- graduated for easy determination of the volume

The additional safety lock prevents uncontrolled opening of the tubes.



free of DNase/RNase,
DNA and ATP

Volume ml	Fig. No.	Packaging	Quantity per Pack	Order No.
0,5	1	2 x 500	1000	56 15 038
1,5	2	1 x 1000	1000	56 15 008
2,0	3	2 x 500	1000	56 15 108

Microcentrifuge Tubes and Micro Centrifuge

free of DNase, RNase
and pyrogens

CE IVD

Microcentrifuge tubes with screw cap and O-Ring Seal

Polypropylene (PP).

- with conical bottom or additional self-standing
- sterilized
- centrifugable up to approx. 25.000 g
- temperature-resistant from -80 °C to +121 °C
- autoclavable up to +121 °C
- Ø x H: 13.1 x 48 mm
- re-sealable bags



Volume ml	Fig. No.	Bottom	Labeling Area	Sterilized	Packaging	Quantity per Pack	Order No.
0.5	1	self-standing	–	■	2 x 500	1000	56 16 021
0.5	2	conical	–	■	2 x 500	1000	56 16 023
1.5	3	self-standing	■	■	2 x 500	1000	56 16 022
1.5	4	conical	■	■	2 x 500	1000	56 16 027
2.0	5	self-standing	■	■	2 x 500	1000	56 16 026
2.0	6	conical	■	■	2 x 500	1000	56 16 025



CE

Micro Centrifuge for PCR Strips and Microcentrifuge Tubes

Space- and energy-saving micro centrifuge with adequate power for routine centrifugations.

- speed 7000/min, max RCF 2680 x g
- lid must be closed for operation mode
- safety: the lid can only be opened when the rotor is not rotating
- rapid start and stop within seconds
- power supply AC100-240V/50Hz/60Hz
- diameter 150 mm, height 117 mm, weight 0.5 kg
- 1 year guarantee
- supplied complete with two rotors:
 - rotor 1 for: 8 microcentrifuge tubes 0.5 ml till 2.0ml
 - rotor 2 for: 16 microcentrifuge tubes 0.2 ml or two pcr strips of 8 tubes 0.2 ml

Product	Quantity per Pack	Order No.
Micro Centrifuge with two rotors	1	71 00 007