



---

Personalized PDF Catalog  
Catalogue Generated January 22, 2019

---

# SCIENCEWARE® 96-Well PCR Tube Rack

## with Cover, Bel-Art

Supplier: Bel-Art Products



Designed for preparing and storing 0.2 mL thin-walled PCR tubes and strip tubes before or after PCR sample processing applications.

- Compact benchtop footprint
- Alphanumeric grid provides easy sample identification
- Color-coded racks allow quick visual identification
- Translucent polypropylene cover

The compact rack features 96 wells that keep tubes securely upright in an 8 x 12 array






Each rack includes translucent polypropylene cover to protect samples and allow viewing without removal. Steam autoclavable at 121°C (250°F).

An alphanumeric grid assists in easy sample identification. Five rack colors assist in sample organization and quick visual identification.

Ordering information: The Assortment pack includes one of each color; flourescent blue, green, pink, yellow, and orange.

### ORDER

VWR Catalog Number	Unit	Price	Quantity	Color	Supplier No.
89230-718	Case of 20	\$414.07		Assortment Pack	F18902-0005
	Pack of 5	\$103.57			
89230-710	Case of 20	\$437.35		Flourescent Blue	F18902-0001
	Pack of 5	\$109.40			
89230-712	Case of 20	\$437.35		Flourescent Green	F18902-0002
	Pack of 5	\$109.40			
89230-720	Case of 20	\$437.35		Flourescent Orange	F18902-0006

	Pack of 5	\$109.40			
89230-714	Case of 20	\$437.35		Flourescent Pink	F18902-0003
	Pack of 5	\$109.40			
89230-716	Case of 20	\$437.35		Flourescent Yellow	F18902-0004
	Pack of 5	\$109.40			

SPECIFICATIONS

Overall Dimensions	12.7 x 8.9 x 3.4 cm (5 x 3 <sup>1</sup> / <sub>2</sub> x 1 <sup>3</sup> / <sub>8</sub> " )
--------------------	--