

NOTE: NALGENE® LDPE, HDPE, PP and PMP containers can be used for long-term storage, but direct UV exposure should be avoided.



A

2250 Autoclavable Carboys with Handles, polypropylene; white polypropylene screw closure, TPE gasket

Popular round shape with handles and 83B closure, autoclavable PP. Excellent for large batches of culture media, distilled water and other solutions. Graduated to contain in 1-gallon and 5-liter increments. Accepts No. 13-1/2 rubber stopper. NOTE: For best results, autoclave using a properly vented closure. See Sterilizing in the Reference section.

Not recommended for use with hazardous materials. Autoclavable/Graduated/Leakproof

Cat. No.2250	-0020	-0050	-0130
Nom. Should. Cap., L	10	20	50
Nom. Should. Cap., gal.	2-1/2	5-1/2	13
Brim Cap., L (approx.)	12.5	24	55
No. per Pkg	1	1	1
No. per Case	6	4	1

*Refer to Bottles/information section for further information on leakproof testing.