

NOTE: NALGENE® LDPE, HDPE, PP and PMP containers can be used for long-term storage, but direct UV exposure should be avoided.



A

2324 Autoclavable Lowboys™, polypropylene; polypropylene spigot and screw closure

Low-profile design easily fits into tight spaces, on shelves or on refrigerators. Can be stacked three high. Good for storing and dispensing stock solutions. Autoclave between uses to eliminate bacterial contamination. Convenient integral handgrip for safe easy transport. Polypropylene offers excellent chemical resistance. Closure size is 63 mm. NOTE: Dismount and disassemble spigot before autoclaving. See Sterilizing in the Reference section. Autoclavable/Leakproof/Patented

| Cat. No.2324 | -0008 | -0015 |
|------------------------|---------------------|---------------------|
| Cap., L | 8 | 15 |
| Cap., gal. | 2 | 4 |
| Approx. H x W x D, mm | 165 x 368 x 305 | 203 x 368 x 368 |
| Approx. H x W x D, in. | 6-1/2 x 14-1/2 x 12 | 8 x 14-1/2 x 14-1/2 |
| No. per Pkg | 1 | 1 |
| No. per Case | 4 | 4 |

Refer to Reference section for further information on leakproof testing



A

2212 Autoclavable Rectangular Carboys, polypropylene; polypropylene screw closure

Space-saving carboys molded in tough, translucent and autoclavable polypropylene. Ideal for storing solutions, and handling large volumes of powders and other solid samples. 3-1/2-inch (8.89-cm) neck opening. Sturdy stainless steel handle attached to molded-in shoulder loops. Calibrated in liters and gallons. NOTE: For best results autoclave using a properly vented closure. Autoclavable/Graduated/Leakproof

| Cat. No.2212 | -0020 | -0050 |
|-----------------|--------------------|---------------------|
| Nom. Cap., L | 9 | 20 |
| Nom. Cap., Gal. | 2 | 5 |
| Dim., mm | 216 x 146 x 356 | 318 x 219 x 381 |
| Dim., in. | 8-1/2 x 5-3/4 x 14 | 12-1/2 x 8-5/8 x 15 |
| No. per Pkg | 1 | 1 |
| No. per Case | 6 | 4 |

Refer to Bottles/Information section for information on leakproof testing



A

2321 Autoclavable Rectangular Carboys with Spigot, polypropylene; polypropylene spigot and screw closure

Autoclavable carboy designed for dispensing solutions, media and sterile water. 3-1/2-inch (8.89-cm) neck opening for easy filling. Sturdy stainless handle. Calibrated to deliver in liters and gallons. Spigot closure (Cat. No. 6423-0010) included. NOTE: Dismount and disassemble spigot before autoclaving. See Sterilizing in the Reference section. Autoclavable/Graduated/Leakproof

| Cat. No.2321 | -0020 | -0050 |
|-----------------|--------------------|---------------------|
| Nom. Cap., L | 9 | 20 |
| Nom. Cap., Gal. | 2 | 5 |
| Dim., mm | 216 x 146 x 356 | 318 x 219 x 381 |
| Dim., in. | 8-1/2 x 5-3/4 x 14 | 12-1/2 x 8-5/8 x 15 |
| No. per Pkg | 1 | 1 |
| No. per Case | 6 | 4 |

Refer to Reference section for further information on leakproof testing