



DESICCATORS "glass" - "non-vacuum"

- Manufactured thick, robust and heat resistant, durable borosilicate glass without vacuum stopcock.
- Excellent smooth grinded body and lid flanges guarantee airtight closure.
- Supplied with a porcelain perforated desiccator plate.



Excellent smooth grinded body and lid flanges guarantee airtight closure.

catalogue number	desiccator diameter	pack quantity
039.01.200	200 mm	1 piece
039.01.250	250 mm	1 piece
039.01.300	300 mm	1 piece



DESICCATORS "glass" - "vacuum"

- Manufactured thick, robust and heat resistant, durable borosilicate glass with a vacuum stopcock.
- Featured with a 2-way high vacuum stopcock enabling consistent and uniform vacuum draw and release.
- Vacuum stopcock can hold maximum vacuum of 1 Torr more than 24 hours.
- Excellent smooth grinded body and lid flanges guarantee airtight closure.
- Supplied with a porcelain perforated desiccator plate.

catalogue number	desiccator diameter	pack quantity
039.02.200	200 mm	1 piece
039.02.250	250 mm	1 piece
039.02.300	300 mm	1 piece



DESICCATORS "P.P & P.C" - "vacuum"

- Manufactured robust polypropylene base with transparent polycarbonate cover with a vacuum stopcock.
- Featured with a 3-way high vacuum stopcock enabling consistent and uniform vacuum draw and release.
- Vacuum stopcock can hold maximum vacuum of 1 Torr more than 24 hours.
- Silicone gasket between body and the lid ensures grease free, airtight closure.
- Supplied with a polypropylene perforated desiccator plate.

catalogue number	desiccator diameter	pack quantity
039.12.200	200 mm	1 piece
039.12.250	250 mm	1 piece
039.12.300	300 mm	1 piece



Various vacuum pumps can be selected from catalogue pages 152.



DESICCATOR LIDS “glass” - “non-vacuum”

- Manufactured thick, robust and heat resistant, durable borosilicate glass without vacuum stopcock.
- Excellent smooth grinded body and lid flanges guarantee airtight closure.


catalogue number	lid diameter	pack quantity
039.06.200	200 mm	1 piece
039.06.250	250 mm	1 piece
039.06.300	300 mm	1 piece



DESICCATOR LIDS “glass” - “vacuum”

- Manufactured thick, robust and heat resistant, durable borosilicate glass with a vacuum stopcock.
- Featured with a 2-way high vacuum stopcock enabling consistent and uniform vacuum draw and release.

catalogue number	lid diameter	pack quantity
039.07.200	200 mm	1 piece
039.07.250	250 mm	1 piece
039.07.300	300 mm	1 piece



DESICCATOR PLATES “porcelain”

- Manufactured from best quality glazed porcelain. They perfectly fit the base of desiccator bodies.
- Featured with perforated holes for perfect heat transfer.

catalogue number	plate diameter	pack quantity
039.05.200	190 mm	1 piece
039.05.250	230 mm	1 piece
039.05.300	300 mm	1 piece

SILICA - GEL DESICCANT BEADS

- Silica-Gel beads are moisture absorbing self - indicating drying agents.
- Perfect substitute of classical desiccator moisture absorbents.
- Beads has a grain size variable from 1 to 3 mm. They are free from any heavy and environmentally compatible.
- Original orange color gradually turns to white as the gel absorbs moisture.
- Color change offers monitoring of the saturation level.
- Silica-Gel beads can be regenerated when heated at a temperature of maximum 120°C.
- The color turns back to original orange color.
- Packed as 1 kg in a sealed P.P bottle in order to prevent preabsorption of moisture.

catalogue number	packing specification
039.16.001	1 kg in P.P bottle

