



**ERLENMEYERS - “clear”
“glass” - “narrow neck”**

- Manufactured from heat resistant, durable borosilicate glass according to ISO 1773 standard.
- Thick wall, heavy base and reinforced rim prevent breakage during washing and transportation.
- Volume graduations and inscriptions are printed with white enamel.

catalogue number	volume	pack quantity
027.01.025	25 ml	10 pieces
027.01.050	50 ml	10 pieces
027.01.100	100 ml	10 pieces
027.01.250	250 ml	10 pieces
027.01.300	300 ml	10 pieces
027.01.500	500 ml	10 pieces
027.01.901	1000 ml	10 pieces
027.01.902	2000 ml	6 pieces
027.01.903	3000 ml	1 piece
027.01.905	5000 ml	1 piece



**ERLENMEYERS - “clear”
“glass” - “screw cap”**

- Manufactured from heat resistant, durable borosilicate glass with a leak-proof screw cap.
- Thick wall, heavy base and reinforced rim prevent breakage during washing and transportation.
- Supplied polypropylene screw cap with large ridges for increased grip.
- Unique internal flanged design of the caps guarantees leak-proof closure.
- Volume graduations and inscriptions are printed with white enamel.

catalogue number	neck size	volume	pack quantity
027.03.100	GL 32	100 ml	1 piece
027.03.250	GL 32	250 ml	1 piece
027.03.500	GL 32	500 ml	1 piece
027.03.901	GL 32	1000 ml	1 piece