



GRADUATED WIDE NECK BOTTLES SERIES "400"

DESCRIPTION

Will withstand considerable pressure before rupture. Secured by tagging points on bottle and cap. Confidential notes can be placed in the space between insert sealing plug and cap. SUITABLE FOR FOODSTUFF ACCORDING TO: • EC Regulation 10/2011

MATERIAL

PE (BODY) PP (CAP)

	Art.	ml	Grad. ml	O. D. Ø mm	h mm	Bocca / Mouth / Col I.D. Ø mm
BOTTIGLIE GRADUATE A COLLO LARGO GRADUATED WIDE NECK BOTTLES FLAcons GRADUÉS À COL LARGE	408	50	10	40	77	24
	409	100	20	48	89	24
	410	250	25	61	126	34,5
	411	500	50	75	158	34,5
	412	1000	100	93	200	45
	413	2000	100	115	247	45



GRADUATED SQUARE BOTTLES WIDE NECK SERIES "600"

DESCRIPTION

Space saving bottles, ideal for long term storage of samples. Approximate graduations allow reproducible filling without the use of cylinders/beakers. Easy grip caps can be wired down for security. SUITABLE FOR FOODSTUFF ACCORDING TO: • EC Regulation 10/2011

MATERIAL

PE (BODY) PP (CAP)

	Art.	ml	Grad.	Dim. mm	Bocca / Mouth / Col I.D. Ø mm
BOTTIGLIE GRADUATE RETTANGOLARI A COLLO LARGO GRADUATED SQUARE BOTTLES WIDE NECK FLAcons GRADUÉS RECTANGULAIRE À COL LARGE	608	25	5	32x 34x 52	18
	609	50	10	38x 38x 70	24
	610	100	20	42x 48x 90	34,5
	611	250	50	57x 60x110	34,5
	612	500	100	70x 80x138	45
	613	1000	100	80x103x176	58
	614	2000	100	100x134x210	58



GRADUATED NARROW NECK BOTTLES SERIES "300"

DESCRIPTION

Double sealing with insert cap to assure tightness. Sloping shoulders reduce tendency for residue deposits. No leaking. Tagging points on shoulder and cap ensures tamper proof closure. Easy grip cap.

SUITABLE FOR FOODSTUFF ACCORDING TO:

- EC Regulation 10/2011

MATERIAL

PE (BODY) PP (CAP)

	Art.	ml	O. D. Ø mm	h mm	Bocca / Mouth / Col I.D. Ø mm	Grad. ml
BOTTIGLIE GRADUATE A COLLO STRETTO GRADUATED NARROW NECK BOTTLES FLACONS GRADUÉS À COL ÉTROIT	317	50	39	76,5	18,5	10
	319	125	46	115	18,5	20
	323	250	60	137	23	25
	324	500	74	165	23	50
	325	1000	93	213	34,5	100
	326	2000	116	273	34,5	100



NARROW NECK BOTTLES - GREY

DESCRIPTION

Completely opaque. Excellent for all photosensitive material whether liquid or solid. Sloping shoulders reduce tendency for residue deposits. No leaking. Insert plug for maximum security. Easy grip cap. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

	Art.	ml	Ø O. D. mm	h mm	Bocca / Mouth / Col I.D. Ø
BOTTIGLIE A COLLO STRETTO - GRIGIO SCURO NARROW NECK BOTTLES - GREY FLACONS À COL ÉTROIT - GRIS FONCÉ	3171	50	39	76,5	18,5
	3191	125	46	115	18,5
	3231	250	60	137	23
	3241	500	74	165	23
	3251	1000	93	213	34,5
	3261	2000	116	273	34,5



NARROW NECK BOTTLES WITH WHITE TAMPER-EVIDENT CAP

DESCRIPTION

The tamper-evident cap incorporates a safety ring which is secured when the cap is screwed on to the bottle. Opening the bottle breaks the safety ring, showing signs of interference. SUITABLE FOR FOODSTUFF.

MATERIAL

HDPE

	Art.	ml	Ø mm	h mm	Bocca / Mouth / Col Ø mm
BOTTIGLIE CON TAPPO A SIGILLO (COLLO Ø 22 MM) NARROW NECK (Ø MM 22) BOTTLES WITH WHITE TAMPER-EVIDENT CAP FLAÇONS À COL ÉTROIT (Ø MM 22) AVEC BOUCHON INVOLABLE BLANC	15700	50	36	84	22
	15701	100	44	101	22
	15702	250	60	136	22
	15703	500	74	170	22
	15704	1000	90	222	22



WIDE NECK BOTTLES WITH BLACK TAMPER-EVIDENT CAP

DESCRIPTION

The tamper-evident cap incorporates a safety ring which is secured when the cap is screwed on to the bottle. Opening the bottle breaks the safety ring, showing signs of interference. SUITABLE FOR FOODSTUFF.

MATERIAL

HDPE

	Art.	ml	Ø mm	h mm	Bocca / Mouth / Col Ø mm
BOTTIGLIE CON TAPPO A SIGILLO (COLLO Ø 30 MM) WIDE NECK (Ø MM 30) BOTTLES WITH BLACK TAMPER-EVIDENT CAP FLAÇONS À COL LARGE (Ø MM 30) AVEC BOUCHON INVOLABLE NOIR	15710	250	60	137	30
	15711	500	74	177	30
	15712	1000	90	230	30
	15713	2000	120	259	45
					Tappo bianco White cap Bouchon blanc



GRADUATED WIDE NECK BOTTLES - PE

DESCRIPTION

Graduated, general purpose bottle meeting Food and Drug regulations with pictograms. Long neck thread ensures leakproof closure. Wide neck for easy filling and emptying of liquid or powder samples. Conforms to DIN 13316 and 168. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

	Art.	ml	Grad.	O. D. Ø mm	h mm	Ø Bocca mm Mouth Ø mm Col I.D. Ø mm	DIN
BOTTIGLIE GRADUATE A COLLO LARGO NORMALIZZATO GRADUATED WIDE NECK BOTTLES FLAcons GRADUÉS À COL LARGE NORMALISÉ	1608	50	10	38	88	24	GL32
	1610	100	20	48	105	24	GL32
	1612	250	25	60	140	38	GL45
	1614	500	100	75	170	38	GL45
	1616	1000	100	95	206	55	GL63
	1618	2000	100	120	252	55	GL63



GRADUATED WIDE NECK BOTTLES - PP

DESCRIPTION

Autoclavable, graduated bottle meeting Food and Drug regulations with pictograms. Very sturdy, ruptureproof bottle. Long neck thread ensures leakproof closure. Wide neck for easy filling and emptying of liquid or powder samples. Conforms to DIN 13316 and 168. SUITABLE FOR FOODSTUFF.

MATERIAL

PP

	Art.	ml	Grad.	O. D. Ø mm	h mm	Ø Bocca mm Mouth Ø mm Col I.D. Ø mm	DIN
BOTTIGLIE GRADUATE A COLLO LARGO NORMALIZZATO GRADUATED WIDE NECK BOTTLES FLAcons GRADUÉS À COL LARGE NORMALISÉ	1621	50	10	38	88	24	GL32
	1623	100	20	48	105	24	GL32
	1625	250	25	60	140	38	GL45
	1627	500	100	75	170	38	GL45
	1629	1000	100	95	206	55	GL63
	1631	2000	100	120	252	55	GL63



GRADUATED NARROW NECK BOTTLES - PP

DESCRIPTION

Autoclavable, graduated, heavy walled bottle. Ideal for transporting toxic waste waters, boiler feed water etc., back to laboratory for analysis as samples can be secured using shoulder and cap tag to avoid sample contamination in transit. Conforms to DIN 13316 and 168. Remove cap before autoclaving. SUITABLE FOR FOODSTUFF.

MATERIAL

PP

	Art.	ml	Grad.	O. D. Ø mm	h mm	Ø Bocca mm Mouth Ø mm Col I.D. Ø mm	DIN
BOTTIGLIE A COLLO STRETTO NORMALIZZATO GRADUATED NARROW NECK BOTTLES FLAcons À COL ÉTROIT NORMALISÉ	1595	50	10	38	92	13	GL18
	1597	100	20	48	108	13	GL18
	1599	250	25	60	150	19	GL25
	1601	500	100	75	182	19	GL25
	1603	1000	100	95	224	23	GL32



GRADUATED NARROW NECK BOTTLES - WITHOUT CAPS

DESCRIPTION

An ideal bottle for water and oils - there are no shoulders to restrict the flow. They are graduated and so can be used as rough measures. Can be secured using shoulder and cap tag to avoid sample contamination in transit. Conforms to DIN 13316 and 168. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

	Art.	ml	Grad.	O. D. Ø mm	h mm	Ø Bocca mm Mouth Ø mm Col I.D. Ø mm	DIN
BOTTIGLIE A COLLO STRETTO NORMALIZZATO, SENZA TAPPO GRADUATED NARROW NECK BOTTLES - WITHOUT CAPS FLAcons À COL ÉTROIT NORMALISÉ, SANS BOUCHON	15801	50	10	38	90	13	GL18
	15821	100	20	48	106	13	GL18
	15841	250	25	60	148	19	GL25
	15861	500	100	75	180	19	GL25
	15881	1000	100	95	222	23	GL32
	15901	2000	500	120	271	23	GL32