

Beakers

Beakers



A

1201 Griffin Low-Form Beakers, polypropylene

For general laboratory use. Superior chemical resistance. All size beakers have easy-to-read silk-screened graduations, plus catalog number, size code, resin code, maximum use temperature and "no flame" symbol. With ribs for easy stacking. Not for use on a hot plate. Meets ISO/DIS 7056 International Standards for Laboratory Plasticware Beakers. Autoclavable/Graduated

Cat. No.1201	-0030	-0050	-0100	-0150	-0250	-0400	-0600	-1000	-2000	-4000
Cap., ml	30	50	100	150	250	400	600	1000	2000	4000
No. per Pkg	12	12	12	12	6	6	4	3	1	1
No. per Case	48	48	48	48	36	36	24	12	6	4



A

1203 Griffin Low-Form Beakers, polymethylpentene

Have the transparency of glass, excellent chemical resistance, and remarkable heat tolerance – can be used continuously at 150°C with intermittent service to 175°C. All size beakers have easy-to-read silkscreened graduations, catalog number, size code, resin code, maximum use to temperature and no flame symbol. Not for use on a hot plate. Meet the ISO/DIS 7056 International Standards for Laboratory Plasticware Beakers. Autoclavable/Transparent/Graduated

Cat. No.1203	-0030	-0050	-0100	-0150	-0250	-0400	-0600	-1000	-2000	-4000
Cap., ml	30	50	100	150	250	400	600	1000	2000	4000
No. per Pkg	12	12	12	12	6	6	4	3	1	1
No. per Case	36	36	36	36	24	24	12	12	4	4



1220 Graduated Beakers with Handle, high-density polyethylene

Molded-in nominal graduations in ounces and milliliters. Rectangular handle for easy carrying and pouring. Extra-heavy wall. Excellent chemical resistance. Not for use on a hot plate. Graduated

Cat. No.1220	-1250	-2500
Cap., ml	1000	3000
Cap., oz.	32	96
No. per Pkg	1	1
No. per Case	6	6



A

1223 Graduated Beakers with Handle, polymethylpentene

Molded-in nominal graduations in ounces and milliliters. Rectangular handle for easy carrying and pouring. Extra-heavy wall. Window-clear PMP. Large rectangular handle. Beakers withstand continuous use at 150°C with intermittent service to 175°C, and are transparent to microwaves. Not for use on a hot plate. Autoclavable/Transparent/Graduated

Cat. No.1223	-1000	-3000
Cap., ml	1000	3000
Cap., oz.	32	96
No. per Pkg	1	1
No. per Case	6	6