



WASH BOTTLE SCREW CAPS

DESCRIPTION

Fit bottles part. no. 15801 - 15821 - 15841 - 15861 - 15881 - 15901.

MATERIAL

PE

TAPPI PER BOTTIGLIE A SPRUZZETTA
WASH BOTTLE SCREW CAPS
BOUCHONS POUR PISSETTES

Art.	DIN STD
3300	GL18
3301	GL25
3302	GL32



WASH BOTTLES

DESCRIPTION

Ideal for distilled water, the stream of liquid is regulated by hand pressure. Dispensing tip gives ultra-fine stream or can be removed to increase flow. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

BOTTIGLIE A SPRUZZETTA
WASH BOTTLES
PISSETTES

Art.	Cap. ml	Ø mm	h mm
179	50	35	96
180	100	43	115
182	250	60	140
183	500	75	180
185	1000	95	220



WASH BOTTLES

DESCRIPTION

The stream of liquid is obtained by blowing through the straight tube. The direction of the stream of liquid can be controlled by hand. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

	Art.	Cap. ml	Ø mm	h mm
BOTTIGLIE A SPRUZZETTA WASH BOTTLES PISSETTES	680	250	60	140
	681	500	75	180
	682	1000	95	220



INTEGRAL WASH BOTTLES

DESCRIPTION

Easy squeeze wash bottles with integral moulded tube. Dispensing tip, which is provided with PE closure cap, can be cut back to increase flow. Cap can be use also with bottles series 15801 - 1631. Special shape for easy grip and stability. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

	Art.	DIN STD	Cap. ml	h mm
BOTTIGLIE A SPRUZZETTA MODELLO INTEGRALE INTEGRAL WASH BOTTLES PISSETTES MODÈLE "INTEGRALE"	1633	GL32	250	140
	1634	GL32	500	178