


DMT-2500 Multi-Tube Vortex Mixer

DMT-2500 Multi-tube vortex mixer can process 50 samples at the same time. Variety of accessories can be selected to meet the mixing needs with different specifications of test tubes.

【 Product Features 】

1. Max can process 50 samples at the same time to make the experiment more convenient and efficient
2. Microprocessor control, simplistic designed appearance, LED display the speed and time.
3. Friendly program designed, built-in short mix and timing operation modes. Operation stable and silent.
4. Large orbital vortex, performance effective. Various types of foam tube frames and pads to option.



Short Mix
Timing Setting
Continuous Work



Time Setting



Most Test tube
Different Kinds Tubes

【 Technical Parameters 】

Model	DMT-2500
Speed Range	500 ~2500rpm
Speed Accuracy	±10rpm
Shaking Orbit	4mm
Time Range	1s ~9999min, short mix, continuous
Max Load	5kgs
Cover Plate Size	184 x 311mm
Voltage	AC 100 ~230V, 50/60Hz
Power	60W
Fuse	250V, 1A, Φ5 x 20
Dimension (WxDxH)	250 x 426 x 480mm
Net Weight	15kgs



【 Optional Accessories 】

Cat#	Description	Dimension
DMT-A	foam tube frame Φ10mm x 50	245 x 132 x 45mm
DMT-B	foam tube frame Φ12mm x 50	245 x 132 x 45mm
DMT-C	foam tube frame Φ13mm x 50	245 x 132 x 45mm
DMT-D	foam tube frame Φ16mm x 50 (suit for 15ml tube)	245 x 132 x 45mm
DMT-E	foam tube frame Φ25mm x 15	245 x 132 x 45mm
DMT-F	foam tube frame Φ29mm x 15 (suit for 50ml tube)	245 x 132 x 45mm
DMT-H	tray pad set (up pad & base pad)	305 x 178.5 x 25mm



DMT-A



DMT-B



DMT-C



DMT-D



DMT-E



DMT-F



DMT-H