

## MS7-S

### 7" Square Magnetic Stirrer

#### Features

- Speed adjustable from 0-1500rpm with a simple dial knob.
- Max. stirring capacity of 10L[H<sub>2</sub>O].
- Ceramic work plate that can prevent chemical erosion with high resistance.



## MS-S

### Magnetic Stirrer

#### Features

- Speed adjustable from 0-1500rpm with a simple dial knob.
- Max. stirring capacity of 20L[H<sub>2</sub>O].
- Stainless steel ceramic work plate which can provide good chemical resistance.



### F100 protective cover

for non-hotplate model,  
use with MS-S



#### Specifications

	MS7-S	MS-S
Work plate dimension[WxD]	184x184mm(7inch)	φ135mm(5 inch)
Work plate material	Glass ceramic	Stainless steel cover with ceramic
Motor type	Shaded pole motor	Brushless DC motor
Motor input power	15W	18W
Motor output power	1.5W	10W
Power	30W	30W
Voltage	100-120/200-240V,50/60Hz	100-240V,50/60Hz
Stirring positions	1	1
Max. stirring quantity[H <sub>2</sub> O]	10L	20L
Max. magnetic bar[length]	80mm	80mm
Speed range	0-1500rpm	0-1500rpm
Speed display	Scale	Scale
Protection class	IP21	IP42
Dimension[WxDxH]	215x360x112mm	160x280x85mm
Weight	3.8kg	2.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH	5-40°C, 80%RH