



---

Personalized PDF Catalog  
Catalogue Generated January 22, 2019

---

# Centrifuge Tube Rack, 15 & 50 mL Globe Scientific


Supplier: Globe Scientific



Rack is designed to hold both 15 and 50mL centrifuge tubes simultaneously

Ships flat and is simple to assemble. The rack accomodates twenty 50mL tubes and thirty 15mL tubes.

## ORDER

VWR Catalog Number	Unit	Price	Quantity	Description	Color	Supplier No.
89498-062	Each	\$21.86		Centrifuge Tube Rack	Orange	456930