

DS2127 Validation Bottle, polycarbonate, white polypropylene screw closure, TPE gasket

Manufactured with all the same materials and molding processes as NALGENE polycarbonate carboys. Allows convenient material compatibility validation for NALGENE polycarbonate carboys, particularly Cat. No. 2251-series and DS2213. Excellent for QA lot sample storage of expensive pharmaceuticals in a container of the same material as the carboys. Autoclavable/Transparent



A

Cat. No.DS2127	-0250
Cap., ml	250
Cap., oz.	8
Closure Size, mm	53B
No. per Case	24

NOTE: NALGENE® LDPE, HDPE, PP and PMP containers can be used for long-term storage, but direct UV exposure should be avoided.

2210 Carboys with Handles, low-density polyethylene; white polypropylene screw closure

A good choice for storage and transport of reagents. Low density polyethylene (LDPE) resin is made using less catalyst than HDPE and is a better choice for reagents sensitive to the polymerization catalyst. Round, graduated to contain in 5L and 1 gallon increments. Wide shoulder handles allow easy carrying and pouring even with gloved hands. Comes with valve-sealing 83B closure and neck will accept no. 13-1/2 rubber stopper. Graduated/Leakproof



Cat. No.2210	-0020	-0040	-0050	-0065	-0130
Nom. Shoulder Cap., L	10	15	20	25	50
Nom. Shoulder Cap., gal.	2-1/2	4	5-1/2	6-1/2	13
Brim Cap., L (approx.)	12.5	18	23	28	54
No. per Case	6	4	4	4	1

2318 Carboys with Spigot, low-density polyethylene; polypropylene spigot and screw closure

Suitable for collecting and dispensing distilled water, reagents, acids. Graduated to deliver in 1-gallon and 5-liter increments. Spigot boss molded in one piece with bottle, for strength and leak resistance. Spigot closure (Cat. No. 6423-0010) included. Four largest sizes (Cat. Nos. 2318-0020, 2318-0050, 2318-0065 and 2318-0130) have convenient wide shoulder handles. Graduated/Leakproof



Cat. No.2318	-0010*	-0020	-0050	-0065	-0130
Nom. Should. Cap., L	4	10	20	25	50
Nom. Should. Cap., gal.	1	2-1/2	5-1/2	6-1/2	13
Brim Cap., L (approx.)	4.1	12	22.5	27.5	54
Closure Size, mm	38-430	83B	83B	83B	83B
No. per Pkg	1	1	1	1	1
No. per Case	6	6	4	1	1

*Cat. No. 2318-0010 is not graduated.