



SHORT SPOUT DROPPING BOTTLE SCREW CAPS

DESCRIPTION

They fit narrow neck standardized bottles part. no. 15801 - 15821 - 15841 - 15861 - 15881 - 15901 and wide neck standardized bottles part. no. 1608 - 1610. Supplied with PE cap for tip closure.

MATERIAL

PE

	Art.	DIN STD	Lungh. beccuccio mm Spout length mm Long. bec mm
TAPPI CON BECCUCCIO CORTO SHORT SPOUT DROPPING BOTTLE SCREW CAPS BOUCHONS À BEC COURT	3303	GL18	16
	3304	GL25	19
	3305	GL32	19



LONG SPOUT DROPPING BOTTLE SCREW CAPS

DESCRIPTION

They fit narrow neck standardized bottles part. no. 15801 - 15821 - 15841 - 15861 - 15881 - 15901 and wide neck standardized bottles part. no. 1608 - 1610. Supplied with PE cap for tip closure.

MATERIAL

PE

	Art.	DIN STD	Lungh. beccuccio mm Spout length mm Long. bec mm
TAPPI CON BECCUCCIO LUNGO LONG SPOUT DROPPING BOTTLE SCREW CAPS BOUCHONS À BEC LONG	3306	GL18	40
	3307	GL25	41
	3308	GL32	42



DROPPING BOTTLES - RANVIER TYPE

DESCRIPTION

Complete with a PE tube with rubber teat. SUITABLE FOR FOODSTUFF.

MATERIAL

PE

Art.	Cap. ml	Ø ext mm	h mm
388	50	35	96
389	100	43	115
391	250	60	140
392	500	75	180