




TEST TUBES - "glass"


➤ Manufactured by automated machines with consistent length and thickness.

➤ **Alternative 2 tubes:** ➤ Made from neutral glass
➤ Made from borosilicate glass

➤ Due to the properties of borosilicate glass highly recommended for over heated applications.

neutral glass test tubes	tube dimensions	wall thickness	pack quantity
077.01.000	6 x 40 mm	0,8/0,9 mm	100 pieces
077.01.001	10x75 mm	0,8/0,9 mm	100 pieces
077.01.021	12x75 mm	0,8/0,9 mm	100 pieces
077.01.002	12x100 mm	0,8/0,9 mm	100 pieces
077.01.003	13x100 mm	0,8/0,9 mm	100 pieces
077.01.031	13x130 mm	0,8/0,9 mm	100 pieces
077.01.004	16x100 mm	0,9/1,0 mm	100 pieces
077.01.005	16x160 mm	0,9/1,0 mm	100 pieces
077.01.006	18x180 mm	0,9/1,0 mm	100 pieces
077.01.007	20x200 mm	0,9/1,0 mm	100 pieces
077.01.008	20x220 mm	0,9/1,0 mm	100 pieces

borosilicate glass test tubes	tube dimensions	wall thickness	pack quantity
077.02.001	10x75 mm	0,9/1,0 mm	100 pieces
077.02.002	12x100 mm	0,9/1,0 mm	100 pieces
077.02.003	13x100 mm	0,9/1,0 mm	100 pieces
077.02.004	16x100 mm	1,0/1,1 mm	100 pieces
077.02.005	16x160 mm	1,0/1,1 mm	100 pieces
077.02.006	18x180 mm	1,0/1,1 mm	100 pieces



TEST TUBES - "glass" - "with screw cap"

➤ Manufactured from neutral glass featured with a polypropylene screw cap.

➤ Hermetic caps guarantee absolute leak-proof closure.

➤ Specially recommended for tissue culture or other similar applications.

➤ Tubes and caps are delivered in separate packings.

catalogue number	tube dimensions	wall thickness	pack quantity
077.03.001	16x100 mm	1,0 / 1,1 mm	100 pieces
077.03.002	16x160 mm	1,0 / 1,1 mm	100 pieces