

Ultrasonic Cleaner DR-DS Series

Derui International Limited				Contact & info.				Tel: 0086-755-33299292					
Derui® Ultrasonic Cleaner								Email: sales@derui-ultrasonic .com					
Customer													
Approval /Cerficate: CE & RoHS & FCC													
Model		Photo						ACCESSORIES					
DR-DS Series													
Feature		1) All the buttons are touch-sensitive 2)Temperature control knob 30-80°C(cleaning liquid temperature setting), Timer control knob 1-30min&4hours(cleaning time setting) 3)START/PAUSE & Memory button, On/Off button. 4) Degas&sweep button *Degassing is the removal of dissolved oxygen trapped in a liquid. In ultrasonic cleaning systems, trapped, dissolved oxygen in the ultrasonic cleaning solution will hinder the ability of the system to clean. This dissolved oxygen must be removed from the system by de-gassing. *Sweep function guarantees an optimum sound field distribution in the cleaning bath due to a modulation of the frequencies.											
Specifications													
Model No.	Ultrasonic power	Heating power	Ultrasonic frequency	Drain	Tank Vol/Size(mm) L×W×D	TransducerPC	Unit Weight	Packing	Packing size(cm)				
DR-DS13	80 W	60 W	37KHz	NO	1.3L/150×138×65	1 x 80w	2.14KG	4 unit/ctn	53x53x35				
DR-DS20	80 W	100 W	37KHz	NO	2L/150×138×100	1 x 80w	2.18KG	4 unit/ctn	53x53x35				
DR-DS22	160 W	100 W	37KHz	NO	2.2L/238×138×65	2 x 80w	2.50KG	4 unit/ctn	49X34X60				
DR-DS30	160 W	150w	37KHz	NO	3L/238×138×100	2 x 80w	3.20KG	4 unit/ctn	49X34X60				
DR-DS40	160 W	220W	37KHz	NO	4L/302×152×100	2 x 80w	3.75KG	2 unit/ctn	43x57x34				
DR-DS45	160 W	220W	37KHz	NO	4.5L/238×138×150	2 x 80w	4.50KG	2 unit/ctn	49x34x39				
DR-DS60	240 W	270W	37KHz	NO	6L/302×152×150	3 x 80w	4.72KG	2 unit/ctn	57x42x39				
DR-DS100	320 W	310W	37KHz	YES	10L/302×239×150	4 x 80w	5.20KG	1 unit/ctn	47x37x43				
DR-DS120	400 W	380W	37KHz	YES	12L/490×138×150	5 x 80w	5.80KG	1 unit/ctn	65x25x47				
DR-DS130	400 W	420W	37KHz	YES	13L/302x239x200	5 x 80w	6.50KG	1 unit/ctn	46x36x47				
DR-DS150	480W	470W	37KHz	YES	15L/325×295×150	6 x 80w	7.50KG	1 unit/ctn	47x42x47				
DR-DS220	480W	470W	37KHz	YES	22L/498×296×150	6 x 80w	8.75KG	1 unit/ctn	65x42x39				
DR-DS280	640 W	620w	37KHz	YES	28L/502×298×200	8 x 80w	10.00KG	1 unit/ctn	65x42x39				