


Water bath shaker



Both spring tray and dedicated clamps available

Description: Thermostatic shaking water bath

Model: SHA-C

Specifications: Cold rolled sheet structure with powder coated, stainless steel cover and interior, Work room: 49*39*16.5cm, universal spring tray size(flask clamps optional with capacity of 12*250ml, or 6*500ml, or 4*1000ml): 41*31cm, Speed: 0-280 times/m, Temperature control range: ambient+5C-100C, Temperature fluctuation: ± 1 C, Heating power: 2000W, Shaking mode: reciprocating(Rotating optional), Amplitude: 20mm, Timer: 0-120m, Digital display for temperature(Digital display for speed optional), with brass drainage faucet

Power Supply: AC 220V/110V 50Hz/60Hz with EU/UK/US/AU/IT/SA Plug

MOQ: 1

Package: unbranded carton (1 unit in 1 box)

G.W/N.W: 40kg/35kg

Package size: 80L×61W×53H (cm)

Special requirements like working room/tray size, clamps capacity/quantity and working mode could be met with custom made models.