



บริการ-ซ่อม
เครื่องมือวิทยาศาสตร์ และเครื่องแก้วทุกชนิด

Lazada Shopee
science.shop

DRAWELL
Artist of Science

LabTech

ANDBON

SCITEK

www.provsci.com

Skyray

MACY

Derui®

Infitek

@provsci



ADD FRIEND

RUWII

SAFOO

YIK
YONGLEKANG



ZHUOYUE



บริษัท โปรวิชั่น ไซแอนติฟิค จำกัด

252 ถนนสวนสยาม แขวงคันนายาว เขตคันนายาว กรุงเทพฯ 10230

โทร : 0 2517 5772-4 / โทรสาร : 0 2517 5772-4 ต่อ 818

252 Suan Siam Road , Khannayao , Khannayao Bangkok 10230

Tel.0 2517 5772-4 / Fax.0 2517 5772-4 - 818 www.provsci.com

All items:Warranty 1 year/Excluded vat 7%
หมายเหตุ : **ทางบริษัทขอสงวนสิทธิ์ ในการเปลี่ยนแปลงราคาและรายละเอียดสินค้าในเอกสารเล่มนี้ (ทุกผลิตภัณฑ์) โดยไม่จำเป็นต้องแจ้งให้ทราบล่วงหน้า**

ANALYTICAL SALE
1 January - 30 June 2024




All items:Warranty 1 year/Excluded vat 7%
หมายเหตุ : **ทางบริษัทขอสงวนสิทธิ์ ในการเปลี่ยนแปลงราคาและรายละเอียดสินค้าในเอกสารเล่มนี้ (ทุกผลิตภัณฑ์) โดยไม่จำเป็นต้องแจ้งให้ทราบล่วงหน้า**

CONTACT



provsci9@gmail.com
 Tel. 0 2517 5772-4
 Provision Scientific

www.provsci.com



SCITEK

ฟรี!

BSC-IIA2-1400

DR-S08

250,000.-

- Biosafety cabinet Class II A2: Negative Air flow 70% Air circulation 30% Exhaust
- 2 filters Supply and Exhaust: AAF ULPA filter
- LCD touch screen
- Internal size 1300x666x600 mm
- External size 1400x836x2100 mm
- Inflow velocity 0.55 m/s



Mini Ductless Fume Hood

Cat.No	FMH-500 LED display	FMH-500II LCD display
External Size (WxDxH)	550x460x700 (mm)	
Internal Size (WxDxH)	480x340x370(mm)	
Airflow Velocity	0.2~0.6 (m/s), Speed adjustable	
UV and LED lamp	33,000.-	35,000.-

Mini Vertical Laminar Flow cabinet

Cat.No	LCB-500- LED display	LCB-500- LCDdisplay
External Size (WxDxH)	550x460x700 mm	
Internal Size (WxDxH)	480x340x370 mm	
Filter	HEPA filter 99.995% at 0.3 micron	
Airflow Velocity	0.2~0.6 (m/s), Speed adjustable	
UV and LED lamp	31,000.-	34,000.-

FMH-800DL

- Air velocity 0.4 – 0.6 m/s
- External size 800x620x1130 mm
- Internal size 780x540x775 mm
- Max opening 650 mm

62,000.-

Mini Vertical Laminar flow

LCB-V800

- External size 710x680x1670 mm
- Internal size 310x590x625 mm
- Air flow 0.2-0.5 m/s with HEPA filter
- UV lamp 8 W at 253.7 nm

63,000.-



ClassIIA2 BSC

BSC-IIA2-2F

- Biosafety cabinet Class II A2: Negative Air flow 70% Air circulation 30% Exhaust
- ขนาด 2 ฟุต
- Inflow velocity 0.52 m/s
- LCD touchscreen.
- Open- Close window with foot switch

125,000.-

UV/Vis SPECTROPHOTOMETER



Cat.No	Optical system	Wavelength	Spectral bandwidth	Wavelength accuracy	SPECIAL PRICE
MC-V1100	Single Beam	352-1000nm	4 nm	± 2.0 nm	25,000.-
MC-UV1700CPC	Split Beam	190-1100nm	1.8 nm	± 0.3 nm	187,000.-
MC-UV1500CSPC	Double Beam	190-1100 nm	2 nm	± 0.5 nm	165,000.-
MC-UV1800CSPC	Double Beam	190-1100 nm	0.5,1,2,4 nm	± 0.3 nm	318,000.-



Ultrasonic Cleaner

CAT.NO	Tank Capacity (Size)	Function	SPECIAL PRICE
DR-DS30	3L / 238x138x100 mm	Degas/Sweep	8,000.-
DR-DS40	4L / 302x152x100 mm	Timer 0-30 min ∞	10,000.-
DR-DS60	6L / 302x152x150 mm	Temp 30-80 °C	15,000.-
DR-DS100	10L / 302x239x150 mm	-	21,000.-
DR-MH30	3L / 238x138x100 mm	Timer 0-30 min ∞	5,000.-
DR-MH40	4L / 150x152x100 mm	Temp 30-80 °C	8,000.-
DR-MH60	6L / 302x152x150 mm	-	11,000.-
DR-MH100	10L / 302x239x150 mm	-	15,000.-
DR-LD30	3L / 238x138x100 mm	Degas	6,500.-
DR-LD40	4L / 302x152x100 mm	LCD Display	9,000.-
DR-LD60	6L / 302x152x150 mm	Timer 0-30 min ∞	13,000.-
DR-LD100	10L / 302x239x150 mm	Temp 30-80 °C	18,000.-



GEM200

- Capacity : 8x200 mm
- Max Sample weight: 3 kg
- Speed : 3000/min at 50 Hz
- Timer : Digital 0 sec – 99 min 59 secs
- Analogue 0-60 Min

1XX,XXX.-

Analogue and Digital Sieve shaker



ตู้เก็บสารเคมี ไร้ท่อ



Acid & alkali storage cabinet

SF-030P

- (PP) Polypropylene
- Double door, Double lock
- (HxWxD) 1120x1090x460
- PP shelf : 1 pcs

28,000.-

Chemical safety cabinet

SF-030

- 114L
- Double door, Ground wire
- (HxWxD) 1120x1090x460
- Steel shelf : 1 pcs

23,700.-

Fire prevention net gas type chemical cabinet

114L

SF-Z030

- Net Gas Type Safety Cabinet
- 3 Point locking
- (HxWxD) 1390x1090x460 mm

ใส่ไส้กรองได้ 8 ชิ้น

105,000.-



Filtered storage cabinet

SF-G819

- Advance voc defector with alarm
- Acrylic glass door
- 798x819x350 mm

ใส่ไส้กรองได้ 2 ชิ้น

66,000.-

Infitek

Vacuum Centrifugal Concentrator

Suitable for concentration of immunoglobulin, concentration of drug metabolites, SPE solid phase extraction, pre-processing and post-processing of liquid chromatography, ADMET/toxicology, polymer chemistry, DNA/RNA purification and concentration, oligomer synthesis, forensic science/drug abuse testing, general laboratory enrichment

- Speed 2000r/min
- Temp.: RT-60 °C
- Final vacuum degree: 0.1 mbar

XXX,XXX.-



LabTech

20 Mbar

DIAPHRAGM VACUUM PUMPS

- Power : 30W
- Speed : 18 L/min
- Vacuum Manual Controller
- Vacuum 20 mbar



33,000.-

VP18R-Plus

Vacuum Pumps



GM-0.33II
Speed 20 L/min

Vacuum
200 mbar

12,000.-



GM-0.5
Positive pressure 30 psi
Speed 30 L/min

Vacuum
200 mbar

13,500.-



GM-0.5II
Speed 30 L/min

Vacuum
50 mbar

16,500.-

Digi Block Digestion System

ED36S

- Max Temperature : 210°C
- Temp accuracy : $\pm 0.2^{\circ}\text{C}$
- Heating material : Aluminum PTFE Coated
- Vessels Material : Customized
- No. Samples : 36X50mL



159,000.-

ED16S

- Max Temperature : 210°C
- Temp accuracy : $\pm 0.2^{\circ}\text{C}$
- Heating material : Aluminum PTFE Coated
- Vessels Material : Customized
- No. Samples : 16X100mL



136,700.-

Separatory funnel shaker



YKD-08

- Vertical shake
- Shake Range: 40 mm Shake Speed: 100-300T/min
- Vessel: 50-2000 mL 8 PCS

128,000.-

YIK
YONGLEKANG

ZYM ultra pure water system (TYPE I WATER)



Capacity 10 L/H Flow rate 1.5-2.0 L/min with pressure tank
Feed water: Tap water, Temp. 5-45°C
RO Water Quality: Conductivity 2-10 µS/cm
Ultrapure Water Quality: Resistivity 18.25 MΩ

1. Microcomputer large animation display LCD screen.
2. 2-way water quality online monitoring system.
3. Time and date settings, timed standby and timed shutdown functions.
4. Excessive emission, various data queries.
5. UP water circulation to extend the life of consumables.

SPECIAL PRICE

Cat.No	Capacity	TOC level	Bacteria	
ZYM-I-10L	Standard	<10 ppb	<1 cfu/mL	11X,XXX.-
ZYM-II-10L	Standard+ Terminal disinfection filter	<10 ppb	<1 cfu/mL	15X,XXX.-
ZYM-III-10L	Standard + Terminal disinfection filter+UV lamp	<5 ppb	<0.1 cfu/mL	16X,XXX.-
ZYM-IV-10L	Standard + Terminal disinfection filter+UV lamp+UF	<5 ppb	<0.1 cfu/mL	19X,XXX.-

ราคานี้ไม่รวมค่าติดตั้ง

Fully Automatic Autoclave



1. Sterilization mode: Liquid, Liquid with warming, Solid, Waste, Agar
2. Sterilization Temp.: 105-138°C Timer: 1-6000 min
3. SMelting Temp.: 60-115°C Timer: 1-6000 min
4. Warming Temp.: 45-79°C Timer: 1-9999 min

Cat.No	Capacity(L)	Chamber Size	Baskets	SPECIAL PRICE
GR60DA	60	400 x 505 mm	2 pcs	175,000.-
GR85DA	85	400 x 700 mm	2 pcs	194,000.-
GR110DA	110	400 x 895 mm	3 pcs	209,100.-

Circulatraler Oil Bath



Cat.No	Temp.(°C)	Tank size(mm)	Flow(L/min)	Pressure(bar)	SPECIAL PRICE
SC-5A	Ambient+8-95	260x140x140	0-15	0.45	19,900.-
SC-15	Ambient+8-100	300x240x200	0-15	0.45	22,400.-
SC-15A	Ambient+8-200	300x240x200	Cycle	-	33,800.-
SC-15B	Ambient+8-200	300x240x200	0-15	0.45	38,300.-

จุดแสดงอุณหภูมิ

จุดแสดงผลความชื้น

เพิ่มความชื้น และเปิดไฟ

ลดความชื้น

ปุ่มตั้งค่า

สัญญาณความชื้น

AD-300S

AD-600S

DS-125S

AD-80S

AD-50S

AD-30S

ANDBON®

ตู้ดูดความชื้นอัตโนมัติ

Dehumidity Range 25%-60% RH
LED Display Automation Control

CAT.NO	CAPACITY	SPECIAL PRICE
AD-30S	30L / 290X320X420mm	4,200.-
AD-50S	50L / 290X320X605mm	5,400.-
AD-80S	80L / 380X390X616mm	7,900.-
DS-125S	125L / 410X390X900mm	13,500.-
AD-300S	300L / 850X430X965mm	21,400.-
AD-600S	600L / 850X1,835X965mm	43,700.-

Sample preparation System

AW-EH-6030

- Program 20 stages of temperature
- Heating material: Graphite coating with Teflon
- Temp. RT-240 °C
- Hole Diameter: 32 mm
- No. of Samples: 60 Samples

AW-P40

- Program 20 stages of temperature
- Heating material: Graphite coating with Teflon
- Temp. RT-240 °C
- Hole Diameter: 16.5 mm

AW-APE-16

- High throughput Vacuum Parallel Evaporator
- No. of Sample: 16
- Sample bottle volume: 250 mL
- Recommend Sample volume: ≤100 mL
- Temp.: RT+5-85 °C

123,000.-

320,000.-

XXX,XXX.-

AWWII